

ActCel T Cell Activation Reagent

Product Name

English Name: ActCel T Cell Activation Reagent

Packaging Specifications

Filling Volume/CatalugueNumber: 1mL/TL-6001-1000

Product Performance

Reactivity Species: Human **Endotoxin:** < 2 EU/mL Appearance: Brown liquid

Intended Use

ActCel T Cell Activation Reagent is synthesized by conjugating anti-human CD3 and CD28 antibodies to a nanomatrix. It enables gentle and efficient activation of T cells or peripheral blood mononuclear cells (PBMCs) for in vitro research applications.

Instructions for Use

Experimental Steps:

1.1 Perform experiments at a 1:200 dilution ratio. Refer to the table below for specific reagent volumes based on culture plate format:

Culture Plate	Maximum Culture Volume	Total T Cell Count	Volume of Reagent per Well
96-well	0.2 mL	0.3 × 10 ⁶	1 μL
48-well	1 mL	1 × 10 ⁶	5 μL
24-well	2 mL	2 × 10 ⁶	10 μL
12-well	4 mL	4 × 10 ⁶	20 μL
6-well	5 mL	5 × 10 ⁶	25 μL

- 1.2 Example for 1×106 T Cells: Resuspend 1×106 T cells in 995 μL X-VIVOTM 15 medium. Add 5 μL ActCel T Cell Activation Reagent (containing 200 IU/mL IL-2). Incubate at 37°C, 5% CO₂ for 48 hours.
- 1.3 After 2-3 days, centrifuge cells at 300g for 10 minutes. Discard supernatant. Add 2 mL fresh complete medium to remove residual reagent. Continue incubation at 37°C, 5% CO₂.
- 1.4 Every 2 days, gently resuspend cells and supplement with fresh complete medium at a 1:2 ratio (e.g., add 4 mL medium to 2 mL original volume on Day 4). Count cells to calculate expansion folds. Maintain culture under standard conditions (37°C, 5% CO₂).

Key Notes:

Website: www.tlbiotechnology.com

- 1. Thoroughly mix the reagent with cells during incubation to enhance activation efficiency.
- 2. Adjust reagent dilution ratio within 1:100–300 based on experimental requirements.



- 3. Centrifugation time for reagent removal may be optimized per user needs.
- 4. Prior to use, centrifuge the product at 1000 rpm for 1 minute to pellet residual beads adhering to the tube cap.

Precautions

This product is intended for *in vitro* cell culture only. Do not use directly in clinical therapy.

Storage Conditions

2-8°C

Expiration Date

6 months

References

Chandler NJ, Call MJ, Call ME. T Cell Activation Machinery: Form and Function in Natural and Engineered Immune Receptors. International Journal of Molecular Sciences. 2020; 21(19):7424. https://doi.org/10.3390/ijms21197424.