

MaxSortin® L-Separation Columns

Product Name

English Name: MaxSortin® L-Separation Columns

Packaging Specifications

Filling Volume/CatalogueNumber: 1 piece / MS-CL01

Product Performance

Total Cell Capacity: 10^7 to 2×10^9 cells

Marked Cell Capacity: 10^5 to 10^8 cells (When separating cells larger than lymphocytes, the capacity of the separation column may decrease)

Void Volume: 1.35 ml

Liquid Storage Volume: 7 ml

Flow Rate: For 1x PBS containing 0.5% BSA (Bovine Serum Albumin): 1.2 to 2.1 ml/min

Applicable Range: Suitable for samples with a diameter of less than 30 μm

Sterility: The separation column is individually packaged aseptically

Intended Use

The MaxSortin® L-Type Separation Column is used in conjunction with the MaxSortin® series of cell sorting magnetic beads and cell sorting buffers. It enables gentle sorting of cells bound by magnetic beads.

Instructions for Use

Experimental Procedure:

- 1.1 Take the MaxSortin® L-Type Separation Column (MS-CL01) and place it on the MACS sorter (130-090-976), and rinse it twice with 1 mL of MaxSortin® cell sorting buffer.
- 1.2 Add the sample, which is obtained after incubating the MaxSortin® series of cell sorting magnetic beads with the cells, into the L-Type Separation Column. Wait for it to flow out naturally. Then, add the cell sorting buffer twice, with 3 mL added each time. Collect the effluent in a 15 mL tube, and the collected cells are the negative selection cells.
- 1.3 After all the cell sorting buffer has flowed out, remove the L-Type Separation Column from the MACS sorter and place it over another new 15 mL centrifuge tube (make a mark). Add 3 mL of the cell sorting buffer to the L-Type Separation Column and use the plunger to directly expel the liquid. The collected cells are the positive selection cells.
- 1.4 Count and perform flow cytometry testing.

Precautions

This product is only applicable to in vitro cell culture and cannot be used directly for clinical treatment.

Storage Conditions

10 to 35°C.

Expiration Date

12 months

Product Use

For research and manufacturing use