Halo Frame

Operating Manual

Why choose Fule?

Our Advantage

The company is a national high-tech enterprise integrating research and development, production and sales of medical devices, with full intelligent processing equipment production line.

Academician and expert studio was established to help Fule improve its R&D capabilities and further deepen production-study-research cooperation; Approved postdoctoral research station.

With complete hardware facilities, excellent research and development team, and close cooperation with clinical experts, we have obtained more than 100 domestic and foreign patents.

Based on the agent cooperation pattern, to establish a nationwide sales and service network, product supply nearly thousand 3 armour hospital in our country, are exported to more than 20 overseas countries.

Content

Product advantage.....01

Instructions for use.....02

Surgical procedure.....03

Product information.....11

Beijing Fule Science&Technology Development Co.,Ltd


Post Code:101204

Sale hot line:010-60999861/2 Fax:010-60999863/8741

Tel:010-60999866/32/75/17 E-mail:YXSC@fulekejicom

Add:No.50, Mafan g West Industry Zone,Pnggu District,Beijing

Surgical procedure

Instructions for Use


- [Indications]**
 - Cobb angle greater than 75°
 - Severe chest wall deformity
 - Severely diminished cardiopulmonary function
 - Severe vertebral rotation, difficult to correct in a single surgery
- [Contraindications]**
 - Spinal stenosis
 - Spondylolisthesis
 - Acute lumbar sprain
 - Tuberculosis, tumors
 - Osteoporosis
 - Lumbar leg pain

Surgical Procedure

[Step 1] Patient Positioning

Install the Halo frame with the patient in the supine position. Choose the appropriately sized head ring, positioned in the "sweatband" area. The entire ring should be 1-2 cm away from the skull at all points.

Halo size


Surgical Procedure

[Step 2] Temporary Fixation

First, perform temporary three-point fixation. Preliminarily adjust the distance and position of the head ring.

Temporary halo fixation


Surgical Procedure

[Step 2] Pin Insertion Points

The anterior pin insertion point on the head ring is located 1 cm above the eyebrow, at the lateral one-third point. If inserting pins under local anesthesia, the patient must keep their eyes closed during placement. If inserted with eyes open, the patient may experience a pulling sensation on the eyelids or even difficulty closing their eyes after the procedure.


Screw placement

Surgical Procedure

[Step 2] Pin Insertion Points


The posterior pin insertion point on the head ring is located 5 cm behind the tip of the ear. The ideal pin placement has the anterior and posterior pins in a diagonal relationship relative to each other.

Surgical Procedure

[Step 2] Tightening

The head ring should be placed horizontally and not tilted. For adults, the head ring requires 4 pins to ensure stability. For children, 6-8 pins are recommended.


Surgical Procedure

[Step 3] Tightening

Tighten all screws, ensuring they are tightened symmetrically. Follow the diagonal sequence shown in the figure: 1, 2, 3, 4. If using a torque wrench, the maximum torque for adults is 0.56kg/cm^3 (Note: Verify unit, likely $\text{kgf}\cdot\text{cm}$ or Nm). For children, reduce by half, using 0.28kg/cm^3 . Be extremely careful not to over-tighten and penetrate the skull.

Halo fixation


Surgical Procedure

[Step 3]

Traction Reduction

If cervical traction reduction is required, the traction ring (or halo) should be positioned at the level of the ears and connected to the head ring. Connect the traction weights to perform gradual closed reduction traction.

Reduction

Surgical Procedure

[Step 4]

Traction Reduction

Position the patient with the head of the bed elevated and the foot lowered. Start with an initial weight of 6 kg, increasing by 1 kg every 6 hours until reaching a maximum weight of 15 kg. It is recommended to take X-rays to check the reduction status after each weight increase.

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